HUMAN CARE

Your Life. Your Way.

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Content

ENGLISH

User Manual	3
Safety Instructions	4
General Description / Intended Use	5
Technical Specifications	6
Product Label	8
Compliances and Standards	9
Suspension Alternatives	9
Hanger Bars for Human Care Lifts	10
Compatibility	11
Sling Application	12
Troubleshooting	14
Care Instructions	15

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User Manual

Important!

You must read the user manual for your device prior to usage. Keep this booklet and information handy for future use.

Symbols

WARNING symbol indicates a potentially hazardous situation which, if not avoided, could result in serious injury to the user or damage to property and/or the device itself.

INFORMATION symbol indicates recommendations and information for proper, trouble-free usage.

WEIGHT CAPACITY symbol indicates maximum user weight for the product. Do not exceed this weight in any circumstance.

CE MARKING symbol indicates product meets requirements of EU Directive 93/42/EEC (MDD) and/or EU Regulation 2017/745 (MDR).

MEDICAL DEVICE symbol indicates that the product is a medical device according to the definition of medical device in EU Directive 93/42/EEC (MDD) and/or EU Regulation 2017/745 (MDR).

RECYCLE according to local regulations.

READ INSTRUCTIONS and ensure you fully understand them before using this product.

Product modification is not permitted.

Before use, always check that the product is not damaged or worn.

Human Care's products are continuously being developed and updated to provide our customers with the highest quality. We reserve the right to make design changes without prior notice.

Always make sure that you have the most recent version of the manual which is available for downloading from our website at www.humancaregroup.com

Contact your local distributor if you have any questions about the product and its use.

Notice to user/patient in case of serious incident


Any serious incident that has occurred in relation to the product, should be reported to the local contact, who reports to the manufacturer, and the competent authority of the country in which the user/patient is established.

Safety Instructions


The process of lifting and moving a person always involves risks. Therefore, you must read the user manuals for the lift and the lifting equipment carefully. Always make sure that the lifting equipment you are using is intended for use with the lift. As a care provider, you are responsible for the client's safety and you must know whether the client will be able to cope with the lifting process. The sling must be fitted to meet the client's specific needs. If you are in any doubt, please contact the supplier.

- Carry out a careful visual inspection of the sling every time you use it.
- Check that the right sling is being used for the client. The model, size and material must meet the client's needs.

Correct

Incorrect

- Plan the moving and lifting process so that it is as safe as possible.
 - Before raising the client, check that the sling is correctly attached to the lift/hanger bar once the lift straps are stretched.
 - It is important that the hooks of the hanger bar are in the correct position, with the opening of the hook facing upwards on both sides to securely hold the straps in place. Recheck the position of the hooks once the lift straps are stretched.
 - It is important to always check that the loops are at the bottom of all hooks (see pictures) and securely placed below the locking pin on all hooks.
- Check that the client is sitting securely in the sling before starting to lift him/her.
 - Check that the lift is balanced when moving, raising and/or lowering the client. If you are using a Human Care lift with telescopic hanger bar, it is important that the bar hooks are in the same position on both sides.
- Our slings have been designed and tested for indoor use with one client at a time and must only be used for transferring a person. This product is not a swing.
 - Make sure no sharp objects are in contact with the sling.

Slings

General Description / Intended Use

Basic Sling

Our Basic sling is a comfortable sling without reinforcements, handles or padding in the back. The sling adapts to the client and offers good support for the entire back. The Basic sling can be used in most lifting situations, fits most clients and is a good choice if it is not established who is going to be using the sling. The Basic sling is offered in thin polyester net (25045) and our unique Soft Air Mesh fabric (25054). Both materials allow moisture and heat to transfer efficiently from the client. This allows the sling to be left in the wheel chair etc. after the transfer is completed. The Soft Air Mesh material is a high friction material that reduces the risk of sliding during lifting and sitting. It also has unique flame-retardant properties. Basic Sling in polyester has an application aid on the lower part of the back. The hand is slid into the open handle, and the thumb is placed on the outside. Then the sling can be more easily controlled and placed correctly on the user when sitting.

If the slings are not required to be left with the client after lifting the Basic sling in our durable polyester fabric (25070) is recommended. The low friction polyester material allows easy application of the sling and is easy to clean. Basic sling has padded leg parts, except from Soft Air Mesh m. 25054. The Basic sling offers an upright sitting position that is suitable for clients with reduced trunk stability but without the need for extra head support.

Low Basic Sling Polyester (25072) and Low Basic Sling Polyester net (25047)

Our Low Basic Sling is a good solution for users who have torso stability and upper body control. The design is similar with the Basic Sling, however it has a shorter back support. The Low Basic Sling provides support from the tailbone up under the arms and under the legs. The arms can be kept inside, but are mostly kept outside of the sling. It is a perfect sling for a more active user.

Low Basic Sling in polyester has an application aid on the lower part of the back. The hand is slid into the open handle, and the thumb is placed on the outside. Then the sling can be more easily controlled and placed correctly on the user when sitting.

The Basic sling is designed for two-, three- and four-point suspension.

Sit Sling and Bath Sling

Our Sit Sling is primarily designed for moving clients between a bed and a wheelchair or for similar purposes. The sling is stable and has a reinforced back support with two removable plastic slats. The client sits in an upright position which makes it easier to transfer him/her to a wheelchair, for example. The sling can be used with the Human Care waist belt 25145/25146. Our Bath Sling is like the Sit Sling. It is made of net and is primarily designed for use in the bath or shower. The back is reinforced with two removable plastic slats. The waist belt cannot be used with this sling.

The Sit Sling and Bath Sling are designed for two-, three- and four-point suspension. The slings can be used with or without the removable plastic slats.

High Basic Sling Polyester (25032) and High Basic Sling Polyester Net (25037)


Our High Basic Sling is designed for situations requiring extra support and can even be used for lifting from a lying position on the floor. The High Basic Sling provides support from the tailbone up over the head and under the legs, which is ideal for users who need a lot of support.

High Basic Sling in polyester has an application aid on the lower part of the back. The hand is slid into the open handle, and the thumb is placed on the outside. Then the sling can be more easily controlled and placed correctly on the user when sitting.

Technical Specifications


It is crucial to use the right size of sling for the client. To help differentiate our sizes, they are colour coded for quick identification. The colour code is found next to the sling's label in the shape of a loop. This loop can also be used to properly store or hang the sling.


Basic Sling

Item	Part No.:	Size	Colour code	A (cm/inch)	B (cm/inch)	C (cm/inch)	Weight Capacity* (kg/lbs)	Weight Capacity (kg/lbs)
Basic Sling Net	25045-JS	Junior small	Purple	22/8,7	62/24,4	50/19,7	15-20/33-44	300/661
Basic Sling Net	25045-JM	Junior medium	Grey	28/11	65/25,6	54/21,3	20-35/44-77	300/661
Basic Sling Net	25045-S	Small	Red	35/13,8	77/30,3	58/22,8	30-60/66-132	300/661
Basic Sling Net	25045-M	Medium	Yellow	40/15,7	88/34,6	63/24,8	55-90/121-198	300/661
Basic Sling Net	25045-L	Large	Green	44/17,3	98/38,6	65/25,6	85-120/187-265	300/661
Basic Sling Net	25045-XL	Extra large	Blue	48/18,9	102/40,2	70/27,6	110-140/243-309	300/661
Basic Sling Net	25045-XXL	Extra extra large	White	52/20,5	108/42,5	74/ 29,1	130-160/287-353	300/661
Basic Sling Soft	25054-JS	Junior small	Purple	22/8,7	62/24,4	50/19,7	15-20/33-44	300/661
Basic Sling Soft	25054-JM	Junior medium	Grey	28/11	65/25,6	54/21,3	20-35/44-77	300/661
Basic Sling Soft	25054-S	Small	Red	35/13,8	77/30,3	58/22,8	30-60/66-132	300/661
Basic Sling Soft	25054-M	Medium	Yellow	40/15,7	88/34,6	63/24,8	55-90/121-198	300/661
Basic Sling Soft	25054-L	Large	Green	44/17,3	98/38,6	65/25,6	85-120/187-265	300/661
Basic Sling Soft	25054-XL	Extra large	Blue	48/18,9	102/40,2	70/27,6	110-140/243-309	300/661
Basic Sling	25070-JS	Junior small	Purple	22/8,7	62/24,4	50/19,7	15-20/33-44	300/661
Basic Sling	25070-JM	Junior medium	Grey	28/11	65/25,6	54/21,3	20-35/44-77	300/661
Basic Sling	25070-S	Small	Red	35/13,8	77/30,3	58/22,8	30-60/66-132	300/661
Basic Sling	25070-M	Medium	Yellow	40/15,7	88/34,6	63/24,8	55-90/121-198	300/661
Basic Sling	25070-L	Large	Green	44/17,3	98/38,6	65/25,6	85-120/187-265	300/661
Basic Sling	25070-XL	Extra large	Blue	48/18,9	102/40,2	70/27,6	110-140/243-309	300/661
Basic Sling	25070-XXL	Extra extra large	White	52/20,5	108/42,5	74/ 29,1	130-160/287-353	300/661
Sit Sling	25030-JS	Junior small	Purple	22/8,7	56/22	53/20,9	15-20/33-44	300/661
Sit Sling	25030-JM	Junior medium	Grey	28/11	63/24,8	58/22,8	20-35/44-77	300/661
Sit Sling	25030-S	Small	Red	35/13,8	76/29,9	62/24,4	30-60/66-132	300/661
Sit Sling	25030-M	Medium	Yellow	40/15,7	86/33,9	62/24,4	55-90/121-198	300/661
Sit Sling	25030-L	Large	Green	44/17,3	98/38,6	63/24,8	85-120/187-265	300/661
Sit Sling	25030-XL	Extra large	Blue	48/18,9	102/40,2	64/25,2	110-140/243-309	300/661
Bath Sling	25035-JS	Junior small	Purple	22/8,7	56/22	53/20,9	15-20/33-44	300/661
Bath Sling	25035-JM	Junior medium	Grey	28/11	63/24,8	58/22,8	20-35/44-77	300/661
Bath Sling	25035-S	Small	Red	35/13,8	76/29,9	62/24,4	30-60/66-132	300/661
Bath Sling	25035-M	Medium	Yellow	40/15,7	86/33,9	62/24,4	55-90/121-198	300/661
Bath Sling	25035-L	Large	Green	44/17,3	98/38,6	63/24,8	85-120/187-265	300/661
Bath Sling	25035-XL	Extra large	Blue	48/18,9	103/3,9	64/25,2	110-140/243-309	300/661


Slings


Low Basic Sling Net

High Basic Sling

Bath Sling

Sit Sling

Item	Part No.:	Size	Colour code	A (cm/inch)	B (cm/inch)	C (cm/inch)	Weight Capacity* (kg/lbs)	Weight Capacity (kg/lbs)
Low Basic Sling Net	25047-JS	Junior small	Purple	22/8,7	65/25,6	30/11,8	15-20/33-44	300/661
Low Basic Sling Net	25047-JM	Junior medium	Grey	28/11	71/28	37/14,6	20-35/44-77	300/661
Low Basic Sling Net	25047-S	Small	Red	35/13,8	83/32,7	42/16,5	30-60/66-132	300/661
Low Basic Sling Net	25047-M	Medium	Yellow	40/15,7	94/37	43/16,9	55-90/121-198	300/661
Low Basic Sling Net	25047-L	Large	Green	44/17,3	104/40,9	45/17,7	85-120/187-265	300/661
Low Basic Sling Net	25047-XL	Extra large	Blue	48/18,9	108/42,5	48/18,9	110-140/243-309	300/661
Low Basic Sling Net	25047-XXL	Extra extra large	White	52/20,5	111/43,7	52/20,5	130-165/287-363	300/661
Low Basic Sling	25072-JS	Junior small	Purple	22/8,7	65/25,6	30/11,8	15-20/33-44	300/661
Low Basic Sling	25072-JM	Junior medium	Grey	28/11	71/28	37/14,6	20-35/44-77	300/661
Low Basic Sling	25072-S	Small	Red	35/13,8	83/32,7	42/16,5	30-60/66-132	300/661
Low Basic Sling	25072-M	Medium	Yellow	40/15,7	94/37	43/16,9	55-90/121-198	300/661
Low Basic Sling	25072-L	Large	Green	44/17,3	104/40,9	45/17,7	85-120/187-265	300/661
Low Basic Sling	25072-XL	Extra large	Blue	48/18,9	108/42,5	48/18,9	110-140/243-309	300/661
Low Basic Sling	25072-XXL	Extra extra large	White	52/20,5	111/43,7	52/20,5	130-165/287-363	300/661
High Basic Sling Net	25037-JS	Junior small	Purple	22/8,7	62/24,4	67/26,4	15-20/33-44	300/661
High Basic Sling Net	25037-JM	Junior medium	Grey	28/11	65/25,6	74/29,1	20-35/44-77	300/661
High Basic Sling Net	25037-S	Small	Red	35/13,8	77/30,3	86/33,9	30-60/66-132	300/661
High Basic Sling Net	25037-M	Medium	Yellow	40/15,7	88/34,6	91/35,8	55-90/121-198	300/661
High Basic Sling Net	25037-L	Large	Green	44/17,3	98/38,6	100/39,4	85-120/187-265	300/661
High Basic Sling Net	25037-XL	Extra large	Blue	48/18,9	102/40,2	102/40,2	110-140/243-309	300/661
High Basic Sling Net	25037-XXL	Extra extra large	White	52/20,5	108/42,5	107/42,1	130-160/287-353	300/661
High Basic Sling	25032-JS	Junior small	Purple	22/8,7	62/24,4	67/26,4	15-20/33-44	300/661
High Basic Sling	25032-JM	Junior medium	Grey	28/11	65/25,6	74/29,1	20-35/44-77	300/661
High Basic Sling	25032-S	Small	Red	35/13,8	77/30,3	86/33,9	30-60/66-132	300/661
High Basic Sling	25032-M	Medium	Yellow	40/15,7	88/34,6	91/35,8	55-90/121-198	300/661
High Basic Sling	25032-L	Large	Green	44/17,3	98/38,6	100/39,4	85-120/187-265	300/661
High Basic Sling	25032-XL	Extra large	Blue	48/18,9	102/40,2	102/40,2	110-140/243-309	300/661
High Basic Sling	25032-XXL	Extra extra large	White	52/20,5	108/42,5	107/42,1	130-160/287-353	300/661

*The weights given are only guidelines and are not a substitute for individual tests. Discrepancies can occur, tolerances ± 1 cm, 0,3".

Product Label

The label on the sling contains the following information:

- Human Care company logo
- This way up arrow and marking the centre line of the sling
- CE and MD marking
- Product name and size
- Picture showing the sling in use
- Part number
- Safe Working Load
- “This side facing out” – outside of the sling, not in touch with the client
- Manufacturer
- Date and country of production
- Instruction labels
 - Do not use fabric softener
 - Visual inspection
 - Latex Free
 - Read the manual before use
- Serial number
- Washing instructions
 - Machine Wash Temperature
 - Do Not Bleach
 - Tumble Dry, Low Heat
 - Do Not Dry Clean
 - Never Iron
- Owner name. (For example, write the client’s name or another unique identifier)
- Indicate bi-annual inspection
- UDI: (01) EAN/GS1/GTIN (11) manufacturing date (21) serial number
- Unique Product Code
- QR code to user manual

HUMAN CARE
Product Group

UP

CE
MD

**Basic Sling,
Medium**

REF 25070-M

SWL < 300kg / < 661lbs

This side facing out

Human Care HC AB
Arstaängsvägen 21B
117 43 Stockholm
Sweden
<http://www.humancaregroup.com>

2023-03-03

SN

335638070016

75°C / 167°F

#

Due date next periodic inspection (every 6 months), mark:

2023	1	2	3	4	5	6	7	8	9	10	11	12
2024	1	2	3	4	5	6	7	8	9	10	11	12
2025	1	2	3	4	5	6	7	8	9	10	11	12
2026	1	2	3	4	5	6	7	8	9	10	11	12

UDI

UPC

(01) 0 7331599 44009 4 (11) 230303 (21) 335638070016

8816086010793

Compliances and Standards

Human Care is an ISO 13485:2016 certified Swedish medical device company. The Quality Management System is in compliance with US 21 CFR part 820.

The product is CE marked in accordance with EU

Directive 93/42/EEC (MDD) and/or EU Regulation 2017/745 (MDR), as class I medical device.

The CE mark is on the product.

This sling has been tested in accordance with ISO 10535.

Suspension Alternatives

Our slings are designed with the premise that all people are different with varying lifting requirements, needs or desires.

An individual assessment of the client's needs must always be carried out before choosing an alternative connection method.

To make all our clients feel as comfortable as possible in our slings, most of them can be used with 2 or 3-point suspensions. Some slings can also be used with a 4-point suspension. Below the options are displayed.

2-point suspension

Suitable for most lifting situations providing a good seating position. 2-point suspension is available with floor lifts, fixed and portable ceiling lifts.

Two-point suspension with crossed leg support: Thread the strap on one leg section through the strap on the other leg section before attaching them to the hanger bar.



Two-point suspension with overlapping leg support: This method is suitable for lifting clients with amputated legs, for example. Extension straps can be attached to the lifting straps to prevent the client from leaning too far back in the seated position.

Two-point suspension with separate leg support: Put the client's legs in each leg section and attach the lifting straps to the hanger bar. Note: This method increases the risk of the client sliding out of the sling

Two-point suspension with Y-connection of the leg parts:

To create, for example, more space around the groin, change the pressure under the user's thighs or adjust the user's leg in relation to the hip, a Y-connection of the leg parts can be used.

Thread the left leg section through the leg support loop that is sewn on the outside of the right leg section, so that the leg parts go parallel to each other. Then attach the left leg part to the left hook of the hanger bar. Then do the same with the right leg part that is attached to the right hook of the hanger bar. This creates a Y connection (the connection looks like a Y from the front, see illustration).

3-point suspension

Offers an elongated sitting position and more even weight distribution. The 3-point suspension is suitable especially for clients carrying an increased risk of pressure sores or who need additional space for the upper body.

Thread the left leg section through the loop on the right leg section. Attach the leg section lifting straps to the leg hooks and the back section lifting straps to the lift.

4-point suspension

A spacious sitting position making it suitable for larger clients or when using a stretcher. 4-point suspension is available with both fixed and portable ceiling lifts or floor lifts.

Hanger Bars for Human Care Lifts

HeliQ
(7275x)

Roomer S / Altair
(55751 45 cm, 17,7" / 55764 55 cm, 21,7")

HeliQ
(92938)

Accessory	Description	Part no.
2-point hanger bar	For stationary lift HeliQ. 3 widths:	35 cm/13,8": 72754 45 cm/17,7": 72755 60 cm/23,6": 72757
2-point hanger bar	For portable lifts Roomer S / Altair, 2 widths:	45 cm/17,7": 55751 55 cm/21,7": 55764
4-point hanger bar	For stationary lift HeliQ	92938
4-point suspension	4-point suspension accessory (49 cm). Enables 4-point lifting with the Altair and Roomer S lifts. The suspension bar is placed on the lift's single leg hook.	55619

Compatibility

Human Care Group slings are designed for use with all Human Care lifts and hanger bars. The products are also compatible with products of other suppliers utilizing the same attachment method for securing slings to a hanger bar and lift. However, due to the wide range of products worldwide, Human Care cannot be responsible for any errors that may result from improper use or application of combinations of slings, hanger bars and lifts of other suppliers. The use of the combined system is at the prescriber's own risk and liability.

To ensure the safety of the client and caregiver, Human Care recommends the following minimum requirements:

1. When using slings with lifts and hanger bars, made by manufacturers other than Human Care, the combined system needs to be designed with the same sling attachment method for the hanger bar of floor lifts and ceiling lifts, meaning loop style slings are combined with a loop style hanger bar, clip style slings are combined with a clip style hanger bar.
2. The involved patient lifts and slings must be CE marked in accordance with Directive 93/42/EEC or Regulation 2017/745 and comply with standard requirements for

patient lifts as stated in EN ISO 10535.

3. The caregiver must follow the manufacturer's instruction manuals and recommendations with respect to the use, care, client, caregiver, cleaning, and inspection of the slings and lifts.
4. Different products on the assembled lift system: i.e., lift, hanger bar, sling, scale and other accessories may have different max. loads. It is always the lowest maximum allowable load indicated for the respective products and the system that applies for the entire system. Always check the markings on all individual products in the system.
5. Individual risk assessment, including a practical test of the system and validation with the specific combination used are required to affirm that the size of the sling is correct for the intended use and compatible to the width and design of the hanger bar. The combined system needs to be documented by the competent assessor.

If there are any questions, please contact your local Human Care office or distributor.

Sling Application

Lifting the client from a lying position

1. Stand to the side of client and turn the client to his/her side. Position the sling so that the centre of the sling lines up with the client's spine. Take note of the marking on the sling "This way out, back". The bottom part of the sling must be directly under the client's coccyx. Lay the client on his/her back. If the bed is height-adjustable, it is a good idea to raise the head of the bed so that the client is in a semi-sitting position. Roll the client over onto his/her other side and pull the sling through. Move the sling down behind the client's back until it reaches the coccyx.

2. Bend the knee and pull the leg sections through between the client's legs. Lower the lift to the height of the client's stomach. When a three-point suspension is used, you should thread the left leg section through the loop on the outside of the right leg section. Attach the leg section lifting straps to the leg hooks on the lift. Attach the back section lifting straps to the lift. Hold the lift away from the client's face with one hand as you raise it above the client's head.


Alternative: The sling can also be used for two- or four-point suspension. In these cases, the loop on the right leg section is not used. Instead the leg sections are crossed.

3. Pause in this position and check that all the lifting straps are securely attached to the hooks.


Slings

4. Continue lifting the client. It may be necessary to lower the bed slightly. Raise the client's legs slightly while turning him/her away from the bed. It is a good idea to put your arm around the client while you are moving him/her to create a greater sense of security.
5. Put the leg sections between the client's knees. Turn the client towards you and fold the sling in, under the client, so that the fabric slides towards you when you pull it out after you have turned back the client.

Lifting the client from a sitting position

1. Pull the sling down behind the client's back. Take note of the marking on the sling "This way out, back". The bottom edge of the sling must be behind the client's coccyx.

2. Pull the leg sections through to the front. Pull the leg sections up between the client's legs. Lower the lift to the height of the client's stomach. When a three-point suspension is used, you should thread the left leg section through the loop on the outside of the right leg section. Attach the leg section lifting straps to the leg hooks on the lift. Attach the back section lifting straps to the lift. Hold the lift away from the client's face with one hand as you raise it above the client's head.

Alternative: The sling can also be used for two- or four-point suspension. In these cases the loop on the right leg section is not used. Instead the leg sections are crossed.

3. Pause in this position and check that all the lifting straps are securely attached to the hooks.

4. Continue lifting the client. It is a good idea to put your arm around the client while you are moving him/her to create a greater sense of security.

5. Fold the leg sections and allow it to slide under itself when you pull out at the outside of the thigh. Pull the sling up behind the client's back.

Troubleshooting

Description	Action
The client is at risk of sliding out	<ul style="list-style-type: none"> The sling may be too large. The sling is not correctly fitted- the length of the loops should be stretched at the same time.
The client is unbalanced when lifted	<ul style="list-style-type: none"> The instructions on the individual test log have not been followed or a new assessment needs to be carried out.
The sling is too small	<ul style="list-style-type: none"> This can be uncomfortable for the client. The sling can cut into the client's flesh and/or his/her back may not have enough support.

Note: It is important that the right sling is used for each client and that it is individually tested to ensure that it meets the client's need and the transfer situation.

Slings

Care Instructions

Check for wear and tear

Check the sling for wear and tear. It is important to check the fabric, the edging, the seams and the lifting straps.

Damaged lifting equipment must not be used.

Cleaning

Human Care slings are intended for individual use. To avoid infections or transmission of bacteria, it is important to clean the slings properly. All slings should be machine washed and dried as soon as possible. Slings should always be cleaned after use in a wet situation.

Always refer to the label of the sling before washing.

Periodic inspection

A periodic inspection must be carried out at least every 6 months. See Human Care's periodic inspection log which can be found on our website at www.humancaregroup.com

Product lifetime

The product's expected lifetime 3 years or 150-180 times washing, if the product is used as intended and maintained according to the manufacturer's instructions, depending on the intensity of use and maximum load applied during use. If the product label is no longer legible, the product should be discarded.

Warranty and support

If you need information or support, please contact your Human Care office, www.humancaregroup.com, or distributor.

Recycling

Once the product has reached the end of its life, it should be disposed according to your country's rules, especially for the separate collection of electrical and electronic products.

Correct disposal helps prevent negative consequences for the environment and human health.

Overview of Sling Material Care		Polyester	Polyester Net	Trevira/Soft Air Mesh	Polyester/Polyester Net	Polyester/Cotton	Ventilated Polyester	Cotton	Polyester Tricot	Coated Nylon
	Machine Wash Maximum Temperature 60°C (140°F) (slings with padding)	●	●							●
	Machine Wash Maximum Temperature 75°C (167°F)	●	●	●	●		●		●	
	Tumble Dry Low, Maximum Temperature 30°C- 54°C (86°F- 130°F)									●
	Tumble Dry Low, Maximum Temperature 60°C (140°F)	●	●	●	●		●		●	
	Do Not Tumble Dry					●		●		
	Do Not Iron	●	●	●	●	●	●	●	●	●
	Do Not Dry Clean	●	●	●	●	●	●	●	●	●
	Do Not Bleach	●	●	●	●	●	●	●	●	●
	Do Not Use Fabric Softener	●	●	●	●	●	●	●	●	●
	Do Not Wash					●		●		
	Dispose of After Use					●		●		
	Disinfection Permitted (see list of approved disinfectants)	●	●		●		●		●	

Disinfection Procedure

Human Care slings are intended to be used with one patient at a time. To avoid possible transmission of infectious diseases, it is important to clean the sling if it has become soiled or if it is to be used by a different patient.

Approved Disinfectants	Solutions to Avoid
Chemical Disinfection: 70% Ethanol, 45% isopropanol or similar	 Steam sterilization
Thermal Disinfection: Wash 65°C (149°F) for a minimum of 10 minutes ¹	 Iodine based disinfectants
Thermal Disinfection: Wash 71°C (160°F) for a minimum of 3 minutes ¹	 Gas sterilization
Note: Add 8 minutes to both thermal disinfection options to ensure heat has penetrated the fabric. Finish with a rinse cycle. ¹	

¹ Ref: HSG (95) 18: Hospital Laundry Arrangements for Used and Infected Linen, Department of Health, NHS (National Health Service, England)

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