



Reliable
Unique Handles
Flexible Frame



The Carl-Oskar is designed for users that find it difficult to maneuver a regular rollator or require extra support. The forearm supports are adjustable in height and angle find the optimal positioning for comfort and stability.



Outstanding Stability

The Carl-Oskar RA is the ideal rollator for users recovering from stroke. It features separate forearm supports and large ball-bearing mounted wheels that ensure walking is easy, even on uneven or rough surfaces. Its reactive frame is durable, even for the more active users. The Carl-Oskar is available with a seat as an option.

grips ensures even distribution of weight. The Carl-Oskar RA can also be equipped with a one-hand brake to allow a user to control hand functions on either the left or right hand side.

Adjustable Forearm Supports

The Carl-Oskar RA features very unique forearm supports that are adjustable in both height and direction. This ensures that the forearm supports can be positioned in the best possible way for user comfort and stability.



Unique Handles

The Carl-Oskar RA features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle tube. The ergonomic design of the

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features



Tray

Practical for resting lightweight items. Maximum load 2kg/4.4lb. Part no. 10601



Splash Guard

Protects legs and clothes. Part no. 10608



Reflector Kit

Includes 30 self-adhesive reflectors. Part no. 91090



Slow-Down Brake

Prevents the rollator from moving too fast and enables a controlled walking speed for the user. Part no.: 10620



Light kit

Illuminate your way while walking. Part no. 91091



Security Lock

Flex cord that protects rollator from theft while unattended. Part no. 91092



Cane holder

Simplifies transport of a cane or umbrella. Part no. 10603




One-Hand Brake

Allows user to control hand functions with one hand (installed left or right). Special order item. Part no. 10653



Curb Climber

Helps raise rollator onto a curb. Part no. 10612



Oxygen Holder

Holds oxygen tube firmly in place. Part no. 10615



Seat Pad

State-of-the-art technology keeps pad cool in summer and warm in winter. Part no. 91093



Basket

Included with rollator in Canada. Ideal for shopping Part no. 10621

Technical Specifications

| | Carl-Oskar RA Low 67 | Carl-Oskar RA Low 72 | Carl-Oskar RA Nor 67 | Carl-Oskar RA Nor 72 |
|------------------------------|--------------------------------|--------------------------------|---------------------------------|---------------------------------|
| Part no. | 12315 | 12415 | 12325 | 12425 |
| Part no. with seat | 12315-1 | 12415-1 | 12325-1 | 12425-1 |
| Handle Height (Min - Max) | 95cm - 113cm/ 37.4" - 44.5" | 95cm - 113cm/ 37.4" - 44.5" | 107cm - 126cm/ 42.1" - 49.6" | 107cm - 126cm/ 42.1" - 49.6" |
| Seat Height | 64cm/25.2" | 64cm/25.2" | 64cm/25.2" | 64cm/25.2" |
| Overall Width | 67cm/26.4" | 72cm/28.3" | 67cm/26.4" | 72cm/28.3" |
| Width between Handles | 47cm/18.5" | 52cm/20.5" | 47cm/18.5" | 52cm/20.5" |
| Distance between Rear Wheels | 52.1cm/20.5" | 52.1cm/20.5" | 52.1cm/20.5" | 52.1cm/20.5" |
| Wheel Size | 30cm x 5cm/ 11.8" x 1.8" | 30cm x 5cm/ 11.8" x 1.8" | 30cm x 5cm/ 11.8" x 1.8" | 30cm x 5cm/ 11.8" x 1.8" |
| Product Weight | 13kg/28.6lb | 13.1kg/28.9lb | 13.2kg/29.1lb | 13.3kg/29.3lb |
| Maximum User Weight | 125kg/275.6lb | 125kg/275.6lb | 125kg/275.6lb | 125kg/275.6lb |

Note: Not all accessories are available for sale in all countries.