

Staffan TW



Table Pad
Indoor / Outdoor
Unique Handles



Highly maneuverable, the Staffan trolley walker is designed for indoor use only. Its narrow overall width makes moving around the home easy and its firm table provides comfort and support for those pain sensitive users.



Outstanding Stability

The Staffan TW is the ideal rollator for pain sensitive users suffering from arthritis, Parkinson's or many other neurological disorders. It features a wide, firm table support and a narrow overall width that makes it easy to move in and around the home.

Fixed Rear Wheels

The Staffan TW features fixed rear wheels and front swivel wheels for enhanced maneuverability as well as a hand operated braking system for enhanced safety. Available in two heights, the Staffan folds easily for storage and travel.



Unique Handles

The Staffan TW features unique front-facing, angle adjustable handles that allow users to find the optimal handle position for comfort and stability. Adjustment is easy with integrated measurements right in the handle tube. The ergonomic design of the grips ensures even distribution of weight. The Staffan TW can also be equipped with a onehand brake to allow a user to control hand functions on either the left or right hand side.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10
info.se@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone: +613.723.6734
Fax: +613.723.1058
info.ca@humancaregroup.com


HUMAN CARE NEDERLAND

Elspeterweg 124
8076 PA Vierhouten
Phone: +31 577 412 171
Fax: +31 577 412 170
info.nl@humancaregroup.com


HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K
Austin, TX 78754
Phone: +512.476.7199
Fax: +512.476.7190
info.us@humancaregroup.com

Accessories/Features




Cane holder
Simplifies transport of a cane or umbrella.
Part no. 10603




Basket
Maximum load 4.5kg/ 10lb. Ideal for shopping.
Part no. 10621



Reflector Kit
Includes 30 self-adhesive reflectors.
Part no. 91090



Oxygen Holder
Holds oxygen tube firmly in place. Part no. 10615



Light kit
Illuminate your way while walking. Part no. 91091



Security Lock
Flex cord that protects rollator from theft while unattended.
Part no. 91092



One-Hand Brake
Allows user to control hand functions with one hand (install left/right). Special order item. Part no. 10653

Technical Specifications

	Staffan TW Low	Staffan TW High
Part no.	43211	43221
Handle Height (Min - Max)	86cm - 108cm / 33.8" - 42.5"	110cm - 125cm / 43.7" - 49.2"
Overall Width	62 cm / 24.4"	62 cm / 24.4"
Wheel Size	20cm x 4cm / 7.9" x 1.6"	20cm x 4cm / 7.9" x 1.6"
Product Weight	11.7kg / 25.8lb	11.9kg / 26.2lb
Maximum User Weight	90 kg / 198.4 lbs	125 kg / 275.6 lb