



HUMAN CARE

Your Life. Your Way.



Bariatric Lifts

Bariatric Lifting Systems

Cost-effective, safe solutions


The EZ bariatric fixed ceiling lift comes with some very unique features, making it the optimal choice for a bariatric fixed lift system. The EZ lift comes with a curtain-pass through which allows the track to bypass existing curtain tracks, saving time and money for installation. With overload protection that ensures a safe lift for those more robust patients, a dual pawl emergency brake system and wide 2" straps, the EZ lift will safely lift and transfer your bariatric patients.

The EZ lift's heavy-duty motor and all-metal drivetrain are designed for extended life with continual usage and its onboard intelligence not only helps monitor usage to help drive consistent safe patient handling initiatives but also captures usage time to help track the frequency of usage. With optional power traverse, continuous charge and inspection reminders, using the lift is as easy as pushing a button.

The EZ lift features an easy-to-read backlit LCD digital display and is proudly manufactured in the United States.



At Human Care we provide solutions to reduce injuries and make safe patient handling possible.



Model Description	Part no.	Dimensions	Product Weight	Weight Capacity
EZ CL1000	72000-NA	36.2cm x 31.1cm x 15.9cm / 14.25" x 12.25" x 6.25"	14kg/31lb	453.4kg/ 1000lb
Controls	Battery	Charger	Motor	
Hand Held remote On unit controls	2-12 V SLA, 5.0Ah	120VAC, 60Hz, 1 Amp Power Cord (NEMA 5)	20 Amp, Electronic soft-start and soft-stop	
Lifting Speed	Lifting Belt	Lift Cycles	Options	
2.3cm/second	2" wide 94" spooled	50 lift cycles per charge at 500 lb load [25 at 1000]	Power traverse, in-rail charging, inspection reminder	



AUSTRALIA

Unit 9, 271 Wells Road
Chelsea Heights, Victoria 3196
Phone:- +61 3 8773 1111
Fax:- +61 3 9773 4008

Email: info.au@humancaregroup.com

CANADA

10-155 Colonnade Road
Ottawa, ON K2E 7K1
Phone:- 613.723.6734
Fax:- 613.723.1058

Email: info.ca@humancaregroup.com

NETHERLANDS

Elspeterweg 124
8076 PA Vierhouten
Phone:- + 31 577 412 171
Fax:- +31 577 412 170

Email: info.nl@humancaregroup.com

UNITED STATES

8006 Cameron Street, Suite K
Austin, TX 78754
Phone:- 512.476.7199
Fax:- 512.476.7190

Email: info.us@humancaregroup.com

SWEDEN

Corporate Office

Årstaängsvägen 21C
117 43 Stockholm
Phone:- +46 8 665 35 00
Fax:- +46 8 665 35 10

Email: info.se@humancaregroup.com



HUMAN CARE

Your Life. Your Way.